

Benchtop High Speed Centrifuge



Description

Benchtop High Speed Centrifuge BHC-1040 small footprint of capacity up to 24x1.5 ml/2.2 ml with speed ranging 16500 rpm. Offer SCT and SBT speed control parameter to avoid mix up and ensure to get maximum samples. Features with LED touch screen with multi-functional knob combination parameters.

Available with 10 profiles for acceleration and for deceleration with high efficiency. For safety we introduced electronic door lock function and, imbalance detection, over-speed protection, inching function, short-spin function accommodates angle rotor to securely hold smaller tube sizes of 12x1.5 ml/2.2 ml and 14000 rpm, 13148xg.

Specifications:

Model	BHC-1040
Shaking Speed	16500 rpm
RCF Speed	18110 × g
Speed Accuracy	±20 rpm
Rotors Capacity	24x1.5 ml/2.2 ml
Materials	Stainless steel centrifugal chamber
Display	LED
Noise	≤60Db
Optional Rotor/ Adaptor	rotors: 16500 rpm, 18110xg, 12x1.5 ml/2.2 ml and 14000 rpm, 16435xg, 48x0.5 ml, 14000 rpm, 13367xg, 24x1.5 ml/2.2 ml
Time Range	1 to 1 hour 39 min/continuation
Power Consumption	80 VA
Dimension	330x390x325 mm
Weight	15 kg
Electrical Requirement	AC 220V/50Hz
Catalogue No.	1401402002

Technical:

Model	BHC-1010	BHC-1020	BHC-1030	BHC-1040
Shaking Speed	16500 rpm			
RCF Speed	18110 × g			
Speed Accuracy	±20 rpm	±10 rpm	±20 rpm	
Rotors Capacity	10x5 ml	8x7 ml	8x10 ml	24x1.5 ml/2.2 ml
Materials	Stainless steel centrifugal chamber			
Display	LED			
Noise	≤60Db			
Optional Rotor/ Adaptor	2 rotors: 16500 rpm, 18110xg, 12x1.5 ml/2.2 ml and 14000 rpm, 13148xg, 10x5 ml	2 rotors: 16500 rpm, 18110xg, 12x1.5 ml/2.2 ml and 14000 rpm, 13148xg, 8x7 ml	2 rotors: 16500 rpm, 18110xg, 12x1.5 ml/2.2 ml and 12000 rpm, 12074xg, 8x10 ml	rotors: 16500 rpm, 18110xg, 12x1.5 ml/2.2 ml and 14000 rpm, 16435xg, 48x0.5 ml, 14000 rpm, 13367xg, 24x1.5 ml/2.2 ml
Time Range	1 to 1 hour 39 min/continuation			
Power Consumption	80 VA			
Dimension	330x390x325 mm			
Weight	15 kg	26 kg	15 kg	
Electrical Requirement	AC 220V/50Hz			

Catalogue No.	1401400502	1401401002	1401401502	1401402002
---------------	------------	------------	------------	------------

Shipping:	2 - 3 weeks Delivery.
speed:	16500 rpm
maximum_temperature:	-
power:	80 VA
capacity:	24x1.5 ml/2.2 ml
model:	BHC-1040

Iris Analytical Ltd.

66 W Flagler St #900
Miami, Florida 33130
United States