

Blood type ID Card Centrifuge DIC-3000



Description

Blood type ID Card Centrifuge DIC-3000 with standard programs are built-in for blood type determination and all kinds experiment of serology, which standardize the test and research of clinical.

Users can set the programs with simple touch screen button. Suitable for blood group testing, micro column gel immune detection.

Spins sample volumes up to 24 cards micro gel rotor with maximum speed 3000 rpm. Driven by DC brushless motor.

Auto-balancing centrifuge for blood bank is particularly designed for the quick centrifuge in blood bank with storages up to 9 programs records for users. For safety we introduced electronic door lock function and, imbalance detection, over-speed protection, inching function, short-spin function.

Specifications:

Model	DIC-3000
Capacity	24 cards
Max RCF	1450xg
Speed Range	3000 rpm
Noise	? 60dB?A?
Timer	1 to 99 min
Weight	16 kg
Electrical Requirement	AC220v/50Hz
Catalogue No.	1402150502

Technical:

Model	DIC-3000
Capacity	24 cards
Max RCF	1450xg
Speed Range	3000 rpm
Noise	? 60dB?A?
Timer	1 to 99 min
Weight	16 kg
Electrical Requirement	AC220v/50Hz
Catalogue No.	1402150502

Shipping:	2 - 3 weeks Delivery.
speed:	3000 rpm
maximum_temperature:	-
power:	-
capacity:	24 cards
model:	DIC-3000

Iris Analytical Ltd.

66 W Flagler St #900
Miami, Florida 33130
United States