

# Cyto Centrifuge CC-1000



## Description

Cyto Centrifuge CC-1000 is an independent device especially designed for cytological applications that transfer cells from fluids to glass slides for microscopic examination.

The CC-1000 has the capacity for up to 0.5 to 4 ml sample holders. Each specimen holder will accommodate a specimen of up to 12 ml. All the cells from a sample are placed on a glass slide evenly in a specific, square area; therefore, screening of the slide is shortened and easy. It is simple to load and unload chambers.

## Specifications:

Model	CC-1000
Capacity	0.5-4 mlx12
Max RCF	730xg
Speed Range	3000 rpm
Dimension	370x432x250 mm
Weight	14 kg
Electrical Requirement	AC220V 50Hz 5A
Catalogue No.	1402300502

## Technical:

Model	CC-1000
Capacity	0.5-4 mlx12
Max RCF	730xg
Speed Range	3000 rpm
Dimension	370x432x250 mm
Weight	14 kg
Electrical Requirement	AC220V 50Hz 5A
Catalogue No.	1402300502

<b>Shipping:</b>	2 - 3 weeks Delivery.
<b>speed:</b>	3000 rpm
<b>maximum_temperature:</b>	-
<b>power:</b>	-
<b>capacity:</b>	0.5-4 mlx12
<b>model:</b>	CC-1000

## **Iris Analytical Ltd.**

66 W Flagler St #900  
Miami, Florida 33130  
United States