

Automatic Decapping Centrifuge



Description

Automatic Decapping Centrifuge ACC-1100 is an innovative and best quality range automatic decapping centrifuge. Suitable for the application in medical and hospitals for blood collection.

Our ACC is best of offering a solution for increased throughput requirements and for integration of new analyzing technologies.

Specifications:

Model	ACC-1100
Capacity	100(75 mm) 84(100 mm)
Max RCF	3030xg
Speed Range	4000 rpm
Noise	? 65dB?A?
Timer	1 s to 99 min, 59 sec
Power Consumption	1100 W
Dimension	470x560x380 mm
Packing Dimension	547x727x422 mm
Weight	50 kg
Electrical Requirement	AC220V/50Hz 6A
Catalogue No.	1402351502

Technical:

Model	ACC-1032	ACC-1048	ACC-1100
Capacity	32(4X8)	48(4X12)	100(75 mm) 84(100 mm)
Max RCF	2980xg	3030xg	3030xg
Speed Range	4000 rpm		
Noise	? 65dB?A?		
Timer	1 s to 99 min, 59 sec		
Power Consumption	370 W	600 W	1100 W
Dimension	470x560x380 mm	420x535x370 mm	470x560x380 mm
Packing Dimension	547x727x422 mm	487x762x422 mm	547x727x422 mm
Weight	38 kg	48 kg	50 kg
Electrical Requirement	AC220V/50Hz 5A	AC220V/50Hz 6A	AC220V/50Hz 6A
Catalogue No.	1402350502	1402351002	1402351502

Shipping:	2 - 3 weeks Delivery.
speed:	4000 rpm
maximum_temperature:	-
power:	1100 W
capacity:	100(75 mm) 84(100 mm)

model:

ACC-1100

Iris Analytical Ltd.

66 W Flagler St #900
Miami, Florida 33130
United States