

Walk-In Fume Hoods FHW-1500



Description

The FHW-1500 Walk-in fume hood supply exhaust system for large laboratory instruments and equipment. Meanwhile, operators can walk into the fume hood to operate the equipment conveniently and directly.

Operators are available to walk into the fume hood to operator if needed.

Large room exhaust system, suitable for large laboratory instruments and equipment.

Technical:

Model		FHW-1000		FHW-1200		FHW-1500	
External Size (W*D*H) mm		1200*850*2200mm		1500*850*2200mm		1800*850*2200mm	
Internal Size (W*D*H) mm		1000*650*1500mm		1300*650*1500mm		1600*650*1500mm	
Internal Size (W*D*H) mm		975		1267		1560	
Max Opening (mm)		1150mm					
Exhaust Airflow Volume		1150m3/h		1500m3/h		1850m3/h	
Inflow Air Velocity		?0.4m/s					
Noise		62dB					
Vibration		XYZ direct ?4?m					
Illumination		?800Lux					
Fluorescent Lamp		21W*2		28W*2		40W*2	
Power Supply		220V±10%, 60/50Hz?110V±10%, 60Hz					
Consumption		400W		500W		500W	
Front Window		Manual, Toughened glass., Anti UV					
Material	Exterior	1.2mm cold-rolled steel with bacteria powder coating.					
	Interior	High grade melamine board with good acid and alkali resistance function					
Standard Accessories		Blower(External) ,Fluorescent lamp, Water tap, Gas tap, PP cup sink, Remote control water valve, Remote control gas valve, Waterproof sockets * 4, ?300mm PVC Duct					
Waterproof Sockets		Total load of 4 waterproof sockets: ?500W					
Package Dimensions (W*D*H) mm		1410x1050x2430mm		1710x1050x2430mm		2010x1050x2430mm	
Exhaust Blower		970*810*580mm		970*810*580mm		970*810*580mm	
Electrical Requirements		110V	220V	110V	220V	110V	220V
Catalogue No.		300250501	300250502	300251001	300251002	300251501	300251502

Shipping:	2 - 3 weeks Delivery.
Fluorescent Lamp:	40W*2
UV Lamp:	-
Max. Opening:	1150mm
power:	500W
capacity:	-

model:	FHW-1500
--------	----------

Iris Analytical Ltd.

66 W Flagler St #900
Miami, Florida 33130
United States