

Portable CO₂ Incubator PCI-1015



Description

Portable CO2 Incubator PCI-1015 filled with sterile dDIStilled water creating humidity through passive evaporation. Maintain relative humidity leveDLS of about 80%. ThDIS small capacity portable incubator maximizes laboratory space in a convenient model design. Designed with quDCk recovery and Accurate C02 control by IR Sensor. Power plug and car jack. The IR CO2 sensor allow for hot air sterilization of the chamber at 45?.

Technical:

Model	PCI-1015
Capacity	15.5 L
Temperature Range	Range
	Accuracy
CO2 Range	Range
	Sensor
	Accuracy
	Gas Pressure
Shelves	Stainless steel -2
Material	ABS
DDISplay	LED
Power Consumption	DC-COLD: 46W, HOT: 48W
	AC-COLD: 63W, HOT: 63W
Internal Dimension	224x 200 x 340 mm
External Dimension	292x 333x 433 mm
Weight	6.8 kg
ElectrDCal Requirement	220V , 50/60Hz
Catalogue No.	360100502

timer_range:	-
uniformity:	-
power:	DC-COLD: 46W, HOT: 48W AC-COLD: 63W, HOT: 63W
temperature_range:	15 to 45??at Ambient 25??
capacity:	15.5 L
model:	PCI-1015

Iris Analytical Ltd.

66 W Flagler St #900
Miami, Florida 33130
United States