

Microplate Reader MR-3600



Description

Microplate Reader MR-3600 has professional software and 3 vibrational linear plate functions for laboratory application. Standard accessories for use with standard 96, and 384 wells plates.

Covers the full wavelength range 400 to 1000 nm and measured as quickly as 6 second per well. Can stored up to 200 groups program and 10000 sample result, more than 10 types of curve fitting equation.

Maintain all important measuring parameters including Endpoint method, kinetic method, inhibition rate, single or dual wavelength measurement function. The built-in multi filter wheel make detection of multiple luminescence wavelengths possible.

Technical:

Model	MR-3200	MR-3400		MR-3600		MR-3800	
Wavelength	400 to 1000 nm	400 to 1000 nm(340 to 1100 nm)				400 to 1000 nm	
Optical Path	8-channels						
Light source	Tungsten Halogen Lamp(22W)						
Optical filter	10 filters(Default configuration 405,450,492,63 0 nm)	12 filters(Default configuration 405,450,492,630 nm)				15 filters(Default configuration 405,450,492,630 nm)	
Accuracy	±0.01 Abs						
Repeatability	0.3%						
Resolution	0.001 Abs						
Display	7-inch LCD	8-inch LCD				Rotating Display 10.4-inch LCD	
Reading Range	0 to 4 Abs					0 to 6 Abs	
Linear Range	0 to 3 Abs					0 to 4 Abs	
Stability	-	±0.003 Abs					
Vibration Plate	-	3 types of linear plate(adjustable 0 to 4 hours)					
Analysis software	Professional software						
Storage	200 groups program and 10000 sample result, more than 10 types of curve fitting equation						
Function	Endpoint method, kinetic method, inhibition rate, single or dual wavelength measurement						
Interface	Wi-Fi						
Input Power							
Dimensions	400x260x200 mm	444x320x318 mm		444x320x190 mm		444x320x220/400 mm	
Weight	6 kg	11 kg		10 kg		11 kg	
Electrical Requirement	DC12V 5A	AC 100/50 Hz	AC 240/60 Hz	AC 100/50 Hz	AC 240/60 Hz	AC 100/50 Hz	AC 240/60 Hz
Catalogue No.	4202005-04	420201001	420201002	420201501	420201502	420202001	420202002

Shipping:	2 - 3 weeks Delivery.
residual_volume:	-
volume:	-

model_type:	96, and 384 wells plates
power:	-
resolution:	0.001 Abs
wavelength:	400 to 1000 nm(340 to 1100 nm)
capacity:	-
model:	MR-3600

Iris Analytical Ltd.

66 W Flagler St #900
Miami, Florida 33130
United States