

Microplate Reader MR-5000



Description

Microplate Reader MR-5000 comes with OSRAM64255 light sources ideal for applications that require accurate positioning, as in the medical and scientific application. An 8-channel manifolds as standard accessories for use with standard 96 wells plates. Covers the absorbance range 0-3.500A and measured quickly, which is significantly faster than other current microplate reader methods. End Point, fixed time, Kinetic make the MR-5000 the best platform for any clinical, biotechnology research and pharmaceutical laboratory applications. Can stored up to 100 program, 1000 patient information, 10000 test records in a single unit. Maintain all important measuring parameters including ABS, Cut-Off, Single Standard, Curve, Multi-Percent, Percent Log, Linear, Exponential, Power, and 4PL Regression. A multi-mode reader has several benefits for laboratories applications ranging from ELISAs to nucleic acid and protein quantitation to cell imaging. Together with multiple microplate services and detection modes into a single, more versatile unit may be an ideal choice especially if you have limited area and performed 12 assays in a single plate.

Technical:

Model	MR-5000	
Absorbance Range	0-3.500A	
Optical Path	0-2.000A	
Measurement Range		
Accuracy	0-2.0A: $\pm 1.0\%$ or $\pm 0.007A$	
Resolution	0-2.0A: $\pm 0.5\%$ or $\pm 0.005A$	
Linearity	$\pm 2.0\%$ or $\pm 0.007A$	
Light Source	OSRAM64255, 8V/20W	
Optical filter	405, 450, 492, 630nm (substitute filters from 400-700nm available on special order)	
Display	5.7" LCD display (320×240 discernibility, 256 gray scale)	
CPU	RISC cpu	
Environment Temperature	0°C-40°C	
Humidity	?85%	
Storage	100 program, 1000 patient information, 10000 test data	
Interface	RS-232C, USB, RJ45, SD interface	
Dimensions	460(×330×190 mm	
Weight	7.5 Kg	
Electrical Requirement	AC110V/50/Hz	220V/60Hz
Catalogue No.	420300501	420300502

Shipping:	2 - 3 weeks Delivery.
residual_volume:	-
volume:	-
model_type:	-
power:	-
resolution:	0-2.0A: 0.5% or 0.005A
wavelength:	-

capacity:	-
model:	MR-5000

Iris Analytical Ltd.

66 W Flagler St #900
Miami, Florida 33130
United States