

Ultra-Low Chest Freezers -164° UCF-5128



Description

Ultra-Low Chest Freezers -164? UCF-5128 build-in foaming technology, and VIP insulation which greatly improves the insulation effect, VIP vacuum insulation board, tightly locks the cold air to ensure the great insulation effect. Features with single-pole lubricated compressor refrigeration technology for better-quality refrigeration capacity. Allow maximum capacity up to 128 L with a temperature range of -110°C to -164? by offering excellent temperature uniformity and stability. Available with upright racks, alarm lamp, voltage compensation, remote alarm, and liquid N2 backup system as optional accessories.

Technical:

Model	UCF-5128
Capacity	128 L
Temperature Range	-110 to -164?
Materials	Internal: structure steel, External: Steel plates with spraying
Alarms	Visual & Audio high and low temperature alarm, High/low voltage, sensor failure alarm, Low battery, Condenser heat dissipation, System failure.
Type	Chest
Evaporator	Coil type evaporator
Display	Digital
Refrigerant	Mixture gas
Casters	6
Compressor	1 pcs
Motor	Fan Motor: EBM
Ambient Temperature	16 to 32°C
Certification	CE
Power Consumption	4600 W
Internal Dimension	490×470×585 mm
External Dimension	1510×985×1120 mm
Net Weight	305 kg
Electrical Requirement	380V/50 Hz
Catalogue No.	5001050503

Shipping:	2 - 3 weeks Delivery.
refrigerant:	Mixture gas
defrosting:	-
model_type:	Chest type
power:	4600W
temperature_range:	-110 to -164?
capacity:	128 L

model:	UCF-5128
--------	----------

Iris Analytical Ltd.

66 W Flagler St #900
Miami, Florida 33130
United States