

# Upright Deep Freezers -40? UDF-260



## Description

Upright Pharmacy Refrigerator 2 to 8? UPR-1130 with transparent door design for easy observation with temperature range 2 to 8?, and forced air refrigeration system with air cooler chiller distribution delivers cold air to different levels within each side of the storage chamber, ensuring consistent cabinet wide cooling with adjustable shelves are adjustable to any height and to fulfill different requests. Ergonomic design offers a clear view of stored items over large framed openings. The thin shape allows for cool recovery of your products, without the concern for swinging door clearances. Customers can choose from two types of interior shelving to suit their needs, security and safety lock to avoid unwanted access, provide storage space for easy sorting of a variety of blood products, and casters design and LED interior light.

## Technical:

Model	UPR-1130	UPR-1236	UPR-1306	UPR-1316	UPR-1406	UPR-1416
Capacity	130 L	236 L	306 L	316 L	406 L	416 L
Type	2-8? Pharmacy Refrigerator					
Cabinet Type	Upright					
Temperature Range	2 to 8 ?					
Cooling	5?					
Alarms	High and low temperature alarms, Sensor error alarm, Power failure alarm, Door ajar alarm					
Refrigerant	R600a	R134a		R600a	R134a	R600a
Cooling Type	Forced Air Cooling					
Defrosting	Auto					
Ambient Temp.			10 to 32?		10 to 32?	
Humidity	<60?					
Certification						
Power Consumption	150 W	150 W	220 W	220 W	230 W	190 W
Internal Dimension	418x390x761 mm	500x420x1090 mm	540x425x1380 mm	540x425x1380 mm	685x425x1380 mm	685x429x1380 mm
External Dimension	500x503x1292 mm	580x560x1785 mm	640x580x1920 mm	640x582x1880 mm	780x580x1920 mm	780x580x1920 mm
Net Weight	53 kg	63 kg	95 kg	92.5 kg	115 kg	120 kg
Gross Weight	57 kg	65 kg	110 kg	100 kg	128 kg	133 kg
Electrical Requirement	220V/50Hz					
Catalogue No.	500350502	500351002	500351502	500352002	500352502	500353002

<b>Shipping:</b>	2 - 3 weeks Delivery.
<b>refrigerant:</b>	R600a
<b>defrosting:</b>	Auto
<b>model_type:</b>	Upright
<b>power:</b>	150 W

<b>temperature_range:</b>	2 to 8 ?
<b>capacity:</b>	130 L
<b>model:</b>	UPR-1130

## Iris Analytical Ltd.

66 W Flagler St #900  
Miami, Florida 33130  
United States