

# Double Pharmacy Refrigerator 2 to 8?



## Description

Double Pharmacy Refrigerator 2 to 8? UPR-2706D Series are 10 shelves refrigerators that are always a spectacular addition to pharmaceutical laboratories, and they are remarkably economical. This allows us a reliable and regular development of our pharmaceutical refrigerators.

Each unit consists of a robust housing and quality constituents provide reliability and durability. Ergonomic design offers a clear view of stored items over large framed openings.

The thin shape allows for cool recovery of your products, without the concern for swinging door clearances. Customers can choose from two types of interior shelving to suit their needs, security and safety lock to avoid unwanted access, provide storage space for easy sorting of a variety of blood products, and castes design and LED interior light.

## Technical:

Model	UPR-2656D	UPR-2706D	
Capacity	656 L	1006 L	
Type	2-8? Pharmacy Refrigerator		
Cabinet Type	Upright		
Temperature Range	2 to 8 ?		
Cooling	5?		
Alarms	High and low temperature alarms, Sensor error alarm, Power failure alarm, Door ajar alarm		
Refrigerant	R600a	CFC-Free	
Cooling Type	Forced Air Cooling		
Defrosting	Auto		
Ambient Temp.	10 to 32?		
Humidity	<60?		
Power Consumption	310 W	472 W	
Internal Dimension	1100x454x1325 mm	1100x684x1325 mm	
External Dimension	1220x642x1885 mm	1220x872x1885 mm	
Net Weight	158 kg	207 kg	
Gross Weight	188 kg	237 kg	
Electrical Requirement	220V/50 Hz	187V/50Hz	242V/60Hz
Catalogue No.	500400502	500401001	500401002

<b>Shipping:</b>	2 - 3 weeks Delivery.
<b>refrigerant:</b>	CFC-Free
<b>defrosting:</b>	Auto
<b>model_type:</b>	Upright
<b>power:</b>	472 W
<b>temperature_range:</b>	2 to 8 ?
<b>capacity:</b>	1006 L
<b>model:</b>	UPR-2706D

**Iris Analytical Ltd.**

**BIOLAB  
SYSTEMS**

66 W Flagler St #900

Miami, Florida 33130

United States