

Low Temperature Chest Freezers -60?



Description

Low-Temperature Chest Freezers -60? LCF-2105C design for large amounts of freezer storage, the available maximum capacities: 105 L with a temperature range of -25°C to -60°C by offering excellent temperature uniformity and stability.

Build-in inner door to avoid more cooling leaking out, and provide more stable inside temperature performance and complete temperature alarm system featuring. Used to stores vaccines, blood plasma, and many other biological products. Features cooling fan for instrument safety, durability, and endurance purpose. The range of medical freezers storage products will permit you to choose the perfect match to meet your requirements through cost and energy saving.

Technical:

| | | |
|------------------------|--|-------------------|
| Model | LCF-2105C | LCF-2485D |
| Capacity | 105 L | 485 L |
| Temperature Range | -25 to 60? | |
| Materials | Inner : Stainless steel; Outer: Color sprayed steel | |
| Alarms | High/low temperature, Sensor error, Thermostat failure | |
| Cooling | Direct cooling | |
| Class | Single door chest | Double door chest |
| Condenser | Copper tube and finned | |
| Sensor | PTC | |
| Refrigerator | CFC-Free | |
| Evaporator | Submerged coil evaporator | |
| Defrosting | Manual | |
| Casters | 4 to 6 | |
| Compressor | CUBIGEL/1 | |
| Temperature Controller | Microprocessor | |
| Ambient Temperature | 10 to 32? | |
| Certification | "ISO9001?13485?14001,CE" | |
| Power Consumption | 280 W | 450 W |
| Internal Dimension | 480x380x665 mm | 1455x525x665 mm |
| External Dimension | 680x580x935 mm | 1655x755x915 mm |
| Net Weight | 62 kg | 124 kg |
| Gross Weight | 70 kg | 138.5 kg |
| Electrical Requirement | 220V/50Hz | |
| Catalogue No. | 500600502 | 500601002 |

| | |
|---------------------------|-----------------------|
| Shipping: | 2 - 3 weeks Delivery. |
| refrigerant: | CFC-Free |
| defrosting: | Manual defrost |
| model_type: | Single door chest |
| power: | - |
| temperature_range: | -25 to -60? |

| | |
|------------------|-----------|
| capacity: | 105 L |
| model: | LCF-2105C |

Iris Analytical Ltd.

66 W Flagler St #900
Miami, Florida 33130
United States