

# Melting Point Apparatus DMP-150



## Description

Digital Melting Point Apparatus DMP-150 convenient and automatic apparatus designed with safety and ease of operation in mind, making it ideal for use in education and industry alike, and touch control display on the front panel incorporates the whole operation. Features with a video camera and a video replay function which allows reviewing the melting point determination. It also suits for color samples and can able to test 4 samples simultaneously, GLP, GMP(EU, UPS & CN), and FDA 21CFR part 11 compliance.

## Technical:

Model	DMP-110	DMP-120	DMP-130	DMP140	DMP-150
Temperature Range	300°C	350°C		400°C	
Position Capacity	4				
Accuracy	±0.4 °C				
Resolution	0.1 °C				
Display	8-Inch touch screen				
Temperature Gradient	Step-less setting at 0.1 to -20°C				
Magnification of Lens	8				
Camera and Video	Camera	Camera And Video			
Video Replay	Manual		Automatic		
Storage	8 GB ROM	16 GB ROM		32 GB ROM	
Graphic	0		100	200	500
Interface	3xUSB, RS232,RJ45, Wi-Fi				
Catalogue No.	5100505-04	5100510-04	5100515-04	5100520-04	5100525-04

<b>Shipping:</b>	2 - 3 weeks Delivery.
<b>interface:</b>	3USB, RS232,RJ45, Wi-Fi
<b>resolution:</b>	0.1 C
<b>temperature_range:</b>	400C
<b>model:</b>	DMP-150

## **Iris Analytical Ltd.**

66 W Flagler St #900  
Miami, Florida 33130  
United States