

Digital Viscometer VM-320



Description

Digital Viscometer VM-320 suitable for the application in pharmaceutical, petroleum, food & beverage industries, pasteurization plant, cosmetic and other personal-care product also industries dealing with paint, rubber, latex, silicon, ink, glue, polymer etc. works on basic principle of fluids internal shear stress and is independent of gravity.

This model come with advanced touch screen with other useful accessories including stepping motor, alarm system, newly design axle with brushless spindles provide smooth operation and easy linear calibration providing high accuracy and precision.

Technical:

Model No.	VM-310		VM-320	
Measuring Range	1 to 100000 mPa.s		1 to 2000000 mPa.s	
Rotor Speed	6, 12, 30, 60 rpm		0.3, 0.6, 1.5, 3, 6, 12, 30, 60 rpm	
Regulator Steps	4		8	
Spindles Type	LV			
Spindle No.	1,2,3,4			
Display	5-inch LCD, Touch screen			
Probe	RTD			
Temperature Range	0 to 100 °C			
Percent Accuracy	±1.0% of range			
Repeatability	±0.5%			
Electrical Requirement	AC100V	AC240V		AC100V AC240V
Catalogue No	650150501	650150502	650151001 650151002	

Shipping:	2 - 3 weeks Delivery.
speed:	0.3, 0.6, 1.5, 3, 6, 12, 30, 60 rpm
reproducibility:	0.5%
measurement_range:	1 to 2000000 mPa.s
temperature_range:	0 to 100 C
model:	VM-320

Iris Analytical Ltd.

66 W Flagler St #900
Miami, Florida 33130
United States