

Portable Viscometer PVM-1010



Description

Portable Digital Viscometer PVM-1010 is highly rugged and portable microcomputer control viscosity meter device used to measure fluid viscosity in laboratory, outdoor and other industrial factory.

Used in the ink, glue, latex, solvent adhesive, polymer solution ,petroleum ,paint ,solvent ,cosmetic ,milk products, pharmaceutical ,juice and so on. Supplied with range of spindles with rotor speed from 60, 100, 150, 200 rpm with measuring range from 200 to 1,200,000 mPa.s with LCD touch screen display.

Technical:

Model No.	PVM-1010	PVM-1020	PVM-1030
Measuring Range	25 to 150,000 mPa.s	50 to 300,000 mPa.s	200 to 1,200,000 mPa.s
Rotor	B1, B2, B3, B4		
Rotor Speed	60, 100, 150, 200 rpm		
Accuracy	±2.0% (full scale)		
Display	LCD		
Repeatability	1.0%		
Electrical Requirement	AC/DC lithium batteries		
Catalogue No	6505005-04	6505010-04	6505015-04

Shipping:	2 - 3 weeks Delivery.
speed:	60, 100, 150, 200 rpm
reproducibility:	1.00%
measurement_range:	25 to 150,000 mPa.s
temperature_range:	-
model:	PVM-1010

Iris Analytical Ltd.

66 W Flagler St #900
Miami, Florida 33130
United States