

Top-loading Balances BT-4500



Description

Top-loading Balances BT-4500 feature weighing capacities of 2100 g, and readabilities of 0.01 g. Some applications include weighing, counting, check weighing, weighing under unstable conditions, peak hold and percentage (%).

Denver top-loading balances are designed for routine weighing applications in your workplace, offering accuracy and repeatability during essential weighing applications in education, laboratory and, industrial sector. It delivers competitive performance at an economical price, featuring RS232 connectivity for easy communication, and a backlit display and a simple interface for uncomplicated operation, these balances are perfectly designed for your workplace.

Standard features include Basic functionality for routine weighing applications. Equipped with three essential weighing modes, RS232 connectivity for data transfer and storage and internal calibration (in selected models) making them ideal for routine weighing applications.

Designed for uncomplicated operation with easy-to-use display and interface Equipped with an easy-to-read, bright backlit display and a simple user interface, they are incredibly easy to operate, with almost no training required.

Smart design and durable construction Denver top pan balances have a small footprint which saves on benchtop space while providing a large weighing surface. They are durably constructed, and feature a stainless-steel pan to withstand day-to-day use in the workplace.

Technical:

Model	BT-4100	BT-4200	BT-4300	BT-4400	BT-4500	BT-4600
Capacity	210 g	310 g	610 g	1100 g	2100 g	3100 g
Repeatability	0.02 g					
Display	LED					
Resolution	0.01 g					
Pan Dimension	130 mm					
Linearity	0.03 g					
Working Temperature	5 to 25°C					
Interface	RS232C, RS485 optional					
Stabilization Time	2 seconds					
Dimension	305x215x115 mm					
Weight	22 kg					
Catalogue No.	68010505-04	68010510-04	68010515-04	68010520-04	68010525-04	68010530-04

Shipping:	2 - 3 weeks Delivery.
pan_dimension:	130 mm
calibration:	-
reproducibility:	0.02 g
resolution:	0.01 g
capacity:	2100 g
model:	BT-4500

Iris Analytical Ltd.

66 W Flagler St #900
Miami, Florida 33130
United States