

Analytical Balances 0.1mg BA-1850



Description

Analytical Balances 0.1mg BA-1850 feature weighing capacities of 310 g, and readabilities of 0.1 mg. Our Analytical Balances are designed with for pharmaceutical, laboratory and regulated weighing in a variety of cell and tissue culture facilities that require the highest accuracy in analytical weighing processes.

Suitable for laboratory, industrial and education application. Particularly it applicable where there are strict hygiene requirements with selectable function including check weighing, parts counting, percentage weighing even the small parts and provide internal calibration. Built-in multiple Units of measurement g/ct/oz make it most convenient for each user. Denver analytical balances have a wide range of readability 0.1 mg, and capacities from 0 to 300 g that cover all of your analytical workflows. Equipped with universal adaptor and round stainless steel weighing pan. All product are CE certified.

Denver Analytical Balances are designed for routine weighing applications in your workplace, offering accuracy and repeatability during essential weighing applications in education, laboratory and, industrial sector. Denver analytical balances deliver competitive performance at an economical price, featuring RS232 connectivity for easy communication, and a backlit display and a simple interface for uncomplicated operation, these balances are perfectly designed for your workplace.

Standard features include Basic functionality for routine weighing applications. Denver analytical balances are equipped with three essential weighing modes, RS232 connectivity for data transfer and storage and internal calibration (in selected models) making them ideal for routine weighing applications.

Designed for uncomplicated operation with easy-to-use display and interface Equipped with an easy-to-read, bright backlit display and a simple user interface, they are incredBTy easy to operate, with almost no training required.

Smart design and durable construction Denver analytical balances have a small footprint which saves on benchtop space while providing a large weighing surface. They are durably constructed, and feature a stainless-steel pan to withstand day-to-day use in the workplace.

Technical:

Model	BA-1050	BA-1150	BA-1250	BA-1350	BA-1450	BA-1550	BA-1650	BA-1750	BA-1850	BA-1950
Capacity	0 to 100 g	0 to 120 g	0 to 160 g	0 to 180 g	0 to 200 g	0 to 210 g	0 to 220 g	0 to 300 g	0 to 310 g	0 to 320 g
Resolution	0.1 mg									
Repeatability	±0.2 mg									
Display	LCD display									
Pan Dimension	80 mm									
Material	Base: Aluminum alloy, Platter: stainless steel									
Linearity	±0.2 mg									
Calibration	Internal calibration									
Interface	RS232									
Certification	CE									
Stabilization Time	3 seconds									
Dimension	340×215×350 mm									
Net Weight	6.8 kg									
Electrical Requirement	NA	AC110V	240V	AC110V	240V	NA				

Catalogue No.	6802005-04	6802010-01	6802010-02	6802015-01	6802015-02	6802020-04	6802025-04	6802030-04	6802035-04	6802040-04	6802045-04	6802050-04
---------------	------------	------------	------------	------------	------------	------------	------------	------------	------------	------------	------------	------------

Shipping:	2 - 3 weeks Delivery.
pan_dimension:	80 mm
calibration:	Internal calibration
reproducibility:	0.2 mg
resolution:	0.1 mg
capacity:	0 to 310 g
model:	BA-1850

Iris Analytical Ltd.

66 W Flagler St #900
Miami, Florida 33130
United States