


Analytical Balances 0.1mg BA-2900



Description

Analytical Balances 0.1mg BA-2900 feature weighing capacities of 310 g, and readabilities of 0.1 mg. Our Analytical Balances are designed with for pharmaceutical, laboratory and regulated weighing in a variety of cell and tissue culture facilities that require the highest accuracy in analytical weighing processes.

Denver Analytical Balances are designed for routine weighing applications in your workplace, offering accuracy and repeatability during essential weighing applications in education, laboratory and, industrial sector. Denver analytical balances deliver competitive performance at an economical price, featuring RS232 connectivity for easy communication, and a backlit display and a simple interface for uncomplicated operation, these balances are perfectly designed for your workplace.

Standard features include Basic functionality for routine weighing applications. Denver analytical balances are equipped with three essential weighing modes, RS232 connectivity for data transfer and storage and internal calibration (in selected models) making them ideal for routine weighing applications.

Designed for uncomplicated operation with easy-to-use display and interface Equipped with an easy-to-read, bright backlit display and a simple user interface, they are incredBTy easy to operate, with almost no training required.

Smart design and durable construction Denver analytical balances have a small footprint which saves on benchtop space while providing a large weighing surface. They are durably constructed, and feature a stainless-steel pan to withstand day-to-day use in the workplace.

Technical:

Model	BA-2100		BA-2200	BA-2300	BA-2400	BA-2500	BA-2600	BA-2700	BA-2800	BA-2900	BA-2980
Capacity	0 to 100 g		0 to 120 g	0 to 160 g	0 to 180 g	0 to 200 g	0 to 210 g	0 to 220 g	0 to 300 g	0 to 310 g	0 to 320 g
Resolution	±0.2 mg										
Display	LCD										
Repeatability	0.1 mg										
Pan Dimension	80 mm										
Material	Base: cast Al Alloy STAINLESS: platter										
Linearity	±0.2 mg										
Calibration	external or internal										
Measuring Parameter	g/mg/ct/oz										
Interface	RS232, USB										
Dimension	310x 212x350 mm										
Weight	5. 0 kg										
Electrical Requirement	AC110V	240V	NA								
Catalogue No.	680250501	680250502	680251004	680251504	680252004	680252504	680253004	680253504	680254004	680254504	680255004

Shipping:	2 - 3 weeks Delivery.
pan_dimension:	80 mm

calibration:	external or internal
reproducibility:	0.1 mg
resolution:	0.2 mg
capacity:	0 to 310 g
model:	BA-2900

Iris Analytical Ltd.

66 W Flagler St #900
Miami, Florida 33130
United States