

Analytical Balances 0.1mg BA-3300



Description

Analytical Balances 0.1mg BA-3300 combines essential weighing functionality with competitive performance, offering high accuracy and reproducibility for applications laboratory, colleges, university, industrial etc. Features with 6 different models with capacities up to 200g with readabilities of 0.1 mg and 80 mm of pan size. Comprise with 5 operational buttons for easy manual operation, automatic internal and external calibration, cast aluminum alloy and stainless steel tray with fast stabilization and response time. Built-in 3 door slide glass offers easy access within the chamber.

Denver Analytical Balances are designed for routine weighing applications in your workplace, offering accuracy and repeatability during essential weighing applications in education, laboratory and, industrial sector. Denver analytical balances deliver competitive performance at an economical price, featuring RS232 connectivity for easy communication, and a backlit display and a simple interface for uncomplicated operation, these balances are perfectly designed for your workplace.

Standard features include Basic functionality for routine weighing applications. Denver analytical balances are equipped with three essential weighing modes, RS232 connectivity for data transfer and storage and internal calibration (in selected models) making them ideal for routine weighing applications.

Designed for uncomplicated operation with easy-to-use display and interface Equipped with an easy-to-read, bright backlit display and a simple user interface, they are incredBTy easy to operate, with almost no training required.

Smart design and durable construction Denver analytical balances have a small footprint which saves on benchtop space while providing a large weighing surface. They are durably constructed, and feature a stainless-steel pan to withstand day-to-day use in the workplace.

Technical:

Model	BA-3100			BA-3200		BA-3300			BA-3400	
Capacity	0 to 100 g			0 to 120 g		0 to 200 g			0 to 220 g	
Resolution	0.1 mg									
Repeatability	±0.2 mg									
Display	LCD									
Pan Dimension	80 mm									
Calibration	External calibration									
Interface	RS232									
Certification	CE									
Stabilization Time	3 to 5 seconds									
Dimension	340×215×350 mm									
Weight	7.5 kg									
Electrical Requirement	AC110V	240V	AC110V		240V		AC110V	240V	AC110V	240V
Catalogue No.	680300501	680300502		680301001	680301002	680301501		680301502	680302001	680302002

Shipping:	2 - 3 weeks Delivery.
pan_dimension:	80 mm
calibration:	External calibration
reproducibility:	0.2 mg

resolution:	0.1 mg
capacity:	0 to 200 g
model:	BA-3300

Iris Analytical Ltd.

66 W Flagler St #900
Miami, Florida 33130
United States