


Touch-Screen Analytical Balances 1mg



Description

Touch-Screen Analytical Balances 1mg BA-5600 feature weighing capacities of 1000 g and readabilities of 1 mg. Our Analytical Balances are designed with for pharmaceutical, laboratory and regulated weighing in a variety of cell and tissue culture facilities that require the highest accuracy in analytical weighing processes.

Suitable for laboratory, industrial and education application. Particularly it applicable where there are strict hygiene requirements with selectable function including check weighing, parts counting, percentage weighing even the small parts and provide automatic internal calibration. Built-in multiple Units of measurement g/mg/ct/oz make it most convenient for each user. Denver premium balances have a wide range of readability 1 mg, and capacities from 0 to 1000 g that cover all of your analytical workflows. Equipped with universal adaptor and round stainless steel weighing pan. Equipped with round stainless steel weighing pan. All product are CE certified. Denver Analytical Balances are designed for routine weighing applications in your workplace, offering accuracy and repeatability during essential weighing applications in education, laboratory and, industrial sector. Denver analytical balances deliver competitive performance at an economical price, featuring RS232 connectivity for easy communication, and a backlit display and a simple interface for uncomplicated operation, these balances are perfectly designed for your workplace.

Standard features include Basic functionality for routine weighing applications. Denver analytical balances are equipped with three essential weighing modes, RS232 connectivity for data transfer and storage and internal calibration (in selected models) making them ideal for routine weighing applications.

Designed for uncomplicated operation with easy-to-use display and interface Equipped with an easy-to-read, bright backlit display and a simple user interface, they are incredBTy easy to operate, with almost no training required.

Smart design and durable construction Denver analytical balances have a small footprint which saves on benchtop space while providing a large weighing surface. They are durably constructed, and feature a stainless-steel pan to withstand day-to-day use in the workplace.

Technical:

Model	BA-5100		BA-5200		BA-5300			BA-5400		BA-5500		BA-5600	
Capacity	0 to 100 g		0 to 200 g		0 to 300g			0 to 500 g		0 to 600 g		0 to 1000 g	
Resolution	1 mg												
Repeatability	±1 mg												
Display	5-inch touch screen												
Pan Dimension	90 mm												
Material	Base: Aluminium alloy, Platter: stainless steel												
Linearity	±1 mg												
Calibration	Automatic Internal calibration												
Interface	RS232												
Certification	CE												
Stabilization Time	3 to 5 seconds												
Dimension	340×215×350 mm												
Weight	7.5 kg												
Electrical Requirement	AC	AC240V	AC	AC240V/DC 5V	AC	AC240V	AC	AC240V	AC	AC240V	AC	AC240V	
	110V/DC 5V	/DC 5V	110V/DC 5V		110V/DC 5V	/DC 5V	110V/DC 5V	/DC 5V	110V/DC 5V	/DC 5V	110V/DC 5V	/DC 5V	

Catalog	680400	6804005	6804010	6804010	680401501	680401	680402	680402	680402	680402	680403	6804030
ue No.	501	02	01	02		502	001	002	501	502	001	02

Shipping:	2 - 3 weeks Delivery.
pan_dimension:	90 mm
calibration:	Automatic Internal calibration
reproducibility:	1 mg
resolution:	1 mg
capacity:	0 to 1000 g
model:	BA-5600

Iris Analytical Ltd.

66 W Flagler St #900
Miami, Florida 33130
United States