

Analytical Balances 1mg BA-7600



Description

Analytical Balances 1mg BA-7600 feature weighing capacities of 320 g, and readabilities of 1 mg. Our Analytical Balances are designed with for pharmaceutical, laboratory and regulated weighing in a variety of cell and tissue culture facilities that require the highest accuracy in analytical weighing processes.

Denver Analytical Balances are designed for routine weighing applications in your workplace, offering accuracy and repeatability during essential weighing applications in education, laboratory and, industrial sector. Denver analytical balances deliver competitive performance at an economical price, featuring RS232 connectivity for easy communication, and a backlit display and a simple interface for uncomplicated operation, these balances are perfectly designed for your workplace.

Standard features include Basic functionality for routine weighing applications. Denver analytical balances are equipped with three essential weighing modes, RS232 connectivity for data transfer and storage and internal calibration (in selected models) making them ideal for routine weighing applications.

Designed for uncomplicated operation with easy-to-use display and interface Equipped with an easy-to-read, bright backlit display and a simple user interface, they are incredBTy easy to operate, with almost no training required.

Smart design and durable construction Denver analytical balances have a small footprint which saves on benchtop space while providing a large weighing surface. They are durably constructed, and feature a stainless-steel pan to withstand day-to-day use in the workplace.

Technical:

Model	BA-7200		BA-7400		BA-7600		BA-7800	
Capacity	120 g		220 g		320 g		420 g	
Repeatability	±0.002 g							
Display	LCD display with backlight and icons							
Resolution	0.001 g							
Pan Dimension	115 mm							
Material	Cast aluminum, static resistant ABS, static-resistant doors, stainless steel pan							
Linearity Error	± 0.002 g							
Calibration	Calibration with External Weight							
Measuring Parameters	g, lb, kg, oz, ct, ozt, pcs, tIs, tIT, tIH, dwt, GN (grain), dr, mom							
Working Temperature	10 °C to 40 °C							
Interface	Configurable RS 232							
Stabilization Time	? 2 sec							
Dimension	304x204x300 mm							
Electrical Requirement	AC110V	AC240V	AC110V	AC240V	AC110V	AC240V	AC110V	AC240V
Catalogue No.	680500501	680500502	680501001	680501002	680501501	680501502	680502001	680502002

Shipping:	2 - 3 weeks Delivery.
pan_dimension:	115 mm
calibration:	Calibration with External Weight
reproducibility:	0.002 g

resolution:	0.001 g
capacity:	320 g
model:	BA-7600

Iris Analytical Ltd.

66 W Flagler St #900
Miami, Florida 33130
United States