

Analytical Balances 1mg BA-8600



Description

Analytical Balances 1mg BA-8600 feature weighing capacities of 200 g, and readabilities of 1 mg. Our Analytical Balances are designed with for pharmaceutical, laboratory and regulated weighing in a variety of cell and tissue culture facilities that require the highest accuracy in analytical weighing processes.

Denver Analytical Balances are designed for routine weighing applications in your workplace, offering accuracy and repeatability during essential weighing applications in education, laboratory and, industrial sector. Denver analytical balances deliver competitive performance at an economical price, featuring RS232 connectivity for easy communication, and a backlit display and a simple interface for uncomplicated operation, these balances are perfectly designed for your workplace.

Standard features include Basic functionality for routine weighing applications. Denver analytical balances are equipped with three essential weighing modes, RS232 connectivity for data transfer and storage and internal calibration (in selected models) making them ideal for routine weighing applications.

Designed for uncomplicated operation with easy-to-use display and interface Equipped with an easy-to-read, bright backlit display and a simple user interface, they are incredBTy easy to operate, with almost no training required.

Smart design and durable construction Denver analytical balances have a small footprint which saves on benchtop space while providing a large weighing surface. They are durably constructed, and feature a stainless-steel pan to withstand day-to-day use in the workplace.

Technical:

Model	BA-8300		BA-8600		BA-8900	
Capacity	0 to 100 g		0 to 200 g		0 to 300 g	
Repeatability	±2 mg					
Display	LED					
Resolution	1 mg					
Pan Dimension	100 mm					
Material	Base: Al alloy, platter: Stainless steel					
Linearity	±2 mg					
Calibration	External					
Stabilization Time	1.5 seconds					
Dimension	290x185x270 mm					
Weight	3.0 kg					
Electrical Requirement	AC110V	AC240V	AC110V	AC240V	AC110V	AC240V
Catalogue No.	680550501	680550502	680551001	680551002	680551501	680551502

Shipping:	2 - 3 weeks Delivery.
pan_dimension:	100 mm
calibration:	External
reproducibility:	2 mg
resolution:	1 mg
capacity:	0 to 200 g
model:	BA-8600

Iris Analytical Ltd.

**BIOLAB
SYSTEMS**

66 W Flagler St #900

Miami, Florida 33130

United States